


The Deep Well Plates from LVL

For the daily use of sample handling and achievement

- Made of clear polypropylene
- SBS standard micro test plate format
- All plates are alphanumerically coded
- Temperature resistant in the range of -80°C to 120°C
- High resistant to solvents
- Plates can be optimally closed with sealing foil or cover mat
- High precision shape, stackable, ideal for processing with automated Pipettor
- Surfaces are suitable for labelling with barcode



1 225.DW.0.3.PP • 2 225.DW.1.2.PP • 3 225.DW.1.2q.PP • 4 225.DW.2.2v.PP • 5 225.DW.5.0.PP • 6 384.DW.0.3.PP

Deep Well Plates

Article No.	Wells	Volume	Cavity	Bottom	Quantity per packaging unit
225.DW.0.3.PP	96	0.3 ml	Round	U-bottom	100
225.DW.1.2.PP	96	1.2 ml	Round	U-bottom	50
225.DW.1.2q.PP	96	1.2 ml	Square	V-bottom	100
225.DW.2.2v.PP	96	2.2 ml	Square	V-bottom	50
225.DW.2.2u.PP	96	2.2 ml	Square	U-bottom	50
225.DW.2.5.PP	96	2.5 ml	Square	U-bottom	20
225.DW.5.0.PP	48	5.0 ml	Rectangular	V-bottom	50
225.DW.10.PP	24	10.0 ml	Rectangular	V-bottom	20
384.DW.0.3.PP	384	0.3 ml	Square	V-bottom	48