

Sterile filter tips

Available volumes:

10 μ l

10 μ l extended

20 μ l

100 μ l

200 μ l

1.000 μ l

While supplies last. Free shipping in Germany on orders over 200€.

- Sterilized
- Certified RNase and DNase free
- Universal fit; compatible with most brands
- Minimum residual wetting
- Easy tip attachment and eject is ideal for single and multichannel pipettes

Order now

1,99€

per Rack

Order now:



Oliver Steinmetz
Sales Manager Germany

oliver.steinmetz@lvl-technologies.com
+49 (0) 172 3525796



www.lvl-technologies.com

Filter Tips for PCR, Molecular Genetics and Biotechnology RNA / DNA FREE

Specially processed filter creates an efficient barrier against aerosols and vapors and protects your pipette from contamination by preventing the penetration of radioactive, environmentally harmful or corrosive substances into the pipette body.

This eliminates the risk of cross-contamination Filter tips have been developed with an optimized cone geometry that ensures a safe and universal fit for all common pipette brands.

	Article numbers	Operating volume	Compatible pipettes						
			Thermo Finpipette®	Brand®	Eppendorf®	Gilson Pipeetman®	Rainin®	DragonLab®	Socorex®
	PX0010-H-ST PXL0010-H-ST	10 µl 10 µl extended	0.2-2, 1-10µL	0.1-2.5, 0.5-10µL	0.1-2.5, 0.5-10µL	0.2-2, 1-10µL	0.5-10µL	0.5-10µL	0.5-10µL
	PX0020-H-ST	20 µl	10-100µL	2-20, 5-50, 10-100µL	2-20µL	2-20µL	2-20µL	2-20µL	2-200µL
	PX0100-H-ST	100 µl	10-100µL	-	10-100µL	10-100µL	10-100µL	10-100µL	20-200µL
	PX0200-H-ST	200 µl	20-200µL	20-200µL	20-200µL	20-200µL	20-200µL	20-200µL	20-200µL
	PX1000-H-ST	1.000 µl	100-1000µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL



Get your free sample box with one rack per variety.

Contact details on the front.



www.lvl-technologies.com