

THE PERFECT SOLUTION FOR SAFE STORAGE AND TRANSPORT OF 100 OR 81 2ML CRYOVIALS!

Practical cryobox for 2ml SAFE® 2D tubes

- **Available in two different sizes:**
100-slot racks for 100 x 2ml cryovials and 81-slot racks for 81 x 2ml cryovials
- **Material:** Durable polycarbonate, temperature-resistant from -198°C to +120°C, autoclavable
- **Lid:** Transparent, non-reversible design, numerically coded
- **Arrangement:** Compartments in 10 x 10 / 9 x 9 rows, numerically labeled
- **Storage:** Suitable for LVL SAFE® 2ml external thread tubes up to -196°C
- **Features:** Bottom holes for plate scanners and tube presenters, linear and HR side barcode, 2D bottom coding can be used as orientation barcode
- **Application:** Ideal for reuse with existing rack towers for LN2 tanks
Perfect for safe storage and easy transport.
Practical cryobox for 2ml SAFE® 2D tubes




SAFE®
Store Aliquots For Ever



Technical Specifications
on the Back





Product Features:

100-Slot / 81-Slot Cryobox

- Material: High quality polycarbonate (PC)
- Temperature resistant: from -198°C to +121°C
- Autoclavable: at 121°C
- Storage: 100x/ 81x 2ml LVL 2D Cryovials

Design and Functionality

- Lid: transparent, made of polycarbonate with non-reversible design
- Bottom: white, with numerical coding
- Arrangement: 10 x 10 compartments/ 9 x 9 compartments, marked numerically
- Fixation: Tubes are firmly fixed, the 2D Floorcode can be scanned
- Area of application: storage at temperatures up to -196°C

Technical Details of the Cryobox:

- **External Dimensions (WxLxH):**
81-Slot Cryo Rack: 132 x 132 x 53,5 mm
100-Slot Cryo Rack: 135 x 135 x 52,5 mm
- Grid Layout: 10 x 10 compartments
- Grid Holes: Optimized for plate scanners and manual tube presenters (TP100/TP81)
- Cryovials: 2ml
- Grid Insert: Permanently integrated
- Material: Polycarbonate (PC)
- **Temperaturbereich:** -196°C bis +121°C
- Features: Autoclavable, sterilizable (Autoclave at +121°C)
- Coding: Numerical coding
- Color: White
- Barcode: Linear and HR side barcode, 2D Bottom coding for orientation

Accessories for cryovials Storage boxes:

For the dry and cool storage, transport, handling, and archiving of samples, we offer you:

- **Stainless steel chest racks**
- **Stainless steel cabinet racks**
- **Stainless steel cabinet inserts**
- **Stainless steel cabinet boxes**

These accessories help improve organization in the laboratory and provide a structured and stable storage solution.



Contact us for more information or visit our website!

SAFE®
Store Aliquots For Ever

LVL
technologies

LVL technologies GmbH & Co. KG
Theodor-Storm-Str. 17 • D-74564 Crailsheim • +49 (0) 7951 956130
info@lvl-technologies.com • www.lvl-technologies.com