


DEFROST THEM ALL!

Now available! LVL SAFE® Advanced Rack Defroster

The device provides a fast and efficient solution to remove ice, making the 2D codes on your sample tubes quickly and reliably readable.



 Technical specifications
located on the back

Find out more:
**LVL Advanced
Rack Defroster**



SAFE 
Store Aliquots For Ever

Order now.

lvl-technologies.com/en/extra/ard/

SAFE® Advanced Rack Defroster

Fast and ergonomic processing of your sample tubes

Our unique electric-powered brush technology quickly de-ices your 2D coded tube rack, ensuring seamless operation in your workflow. Simply press the activation button and slide the frosted rack along the track where the precisely positioned brush removes ice efficiently. This process clears the tube bottoms of any frost or ice residues, enabling flawless and error-free 2D code reading.

No need for chemicals or manual brushing—this innovative solution enhances sample processing efficiency, convenience and reliability in any laboratory setting.

- **Fast & Efficient:** Removes ice in no time
- **Chemical-Free:** Sustainable de-icing solution
- **Technical:** Safe for plastic materials
- **Ergonomic Design:** Comfortable & easy to use

Technical details

Article number	ARD-SBS
Dimensions (L x W x H)	266 x 200 x 86 mm
Weight	4785 g
Power	110 V / 230 V / 50 Hz
Compatible	with all 2D SAFE® tube racks
Standardized	SBS format