


# READ THEM ALL!

## Now available! LVL SAFE® Cryo Reader

Reads cryoboxes 9 x 9, 10 x 10, 14 x 14 and all commonly used SBS racks in less than a second!

- Scan all boxes of racks with one reader.
- No glass: no fog, no reflection, no smudges.
- Power supply: No additional cable required - connect your scanner via USB cable and off you go.



Technical specifications  
located on the back

**SAFE** ❄️  
Store Aliquots For Ever

# Order now.


[www.lvl-technologies.com](http://www.lvl-technologies.com)

## SAFE® Cryo Reader

### Specs & features:

- **Robot/Integration Friendly:** Easily set up the Reader with robotics and automation systems.
- **Adaptability:** Three top plates included for 9x9, 10x10, 14x14 & SBS racks.
- **No Mode Changes Required:** Switch between different racks and tubes instantaneously without the need for reconfiguration - software adjusts to optimal reading parameters.
- **LIMS Integration:** Compatible with most available LIMS systems.
- **GxP ready:** Including User management and audit trail
- **No license key required:** Can be connected to an unlimited number of different PCs.

for:



### Technical details

Article number	DMTR-CRYO
Dimensions (D x W x H)	144 x 163 x 177 mm
Weight	2,25 kg
Supported Operating System	Windows 10 & 11
Safety Approvals	Conforms CE & FCC