

READ THEM ALL!

Now available! LVL SAFE® Integration Reader

Scan SBS Format Sample Racks - decodes any 24, 48, 96, 138, 384 in SBS format well racks without any configuration in less than a second!

- The Integration Reader is engineered for automated lab systems, ensuring precise fit, high compatibility, and streamlined integration.
- No glass: no fog, no reflection, no smudges.
- **Power supply:** No additional cable required - connect your scanner via USB cable and off you go.



Technical specifications
located on the back

Order now.

www.lvl-technologies.com

SAFE® Integration Reader

Specs & features:

- **Robot/Integration Friendly:** Easily set up the Reader with robotics and automation systems.
- **No Mode Changes Required:** Switch between different racks and tubes instantaneously without the need for reconfiguration - software adjusts to optimal reading parameters.
- **LIMS Integration:** Compatible with most available LIMS systems.
- **GxP ready:** Including User management and audit trail
- **No license key required:** Can be connected to an unlimited number of different PCs.
- **Multiple Export Options:** CSV, JSON, PDF, XLS, XMS & TCP/IP data export options available.

for example usable for:



Technical details

Article number	DMTR-INTE
Dimensions (D x W x H)	85 x 150 x 135 mm
Weight	1 kg
Supported Operating System	Windows 10 & 11
Safety Approvals	Conforms CE & FCC